

Strength of Vinegar by Acid-Base Titration

Test Exercise - 100 points

1

? QUESTIONS ?

How are acid/base titrations conducted?

What is standardization?

How do you standardize a solution of a base?

How, and with what accuracy and precision can you determine the concentration of a solution of NaOH?

Using that NaOH solution, with what accuracy and precision can you determine the concentration of a solution of acetic acid?

2

Concepts:

Strong/Weak Acids Acid Dissociation/ K_a End point
Equivalence point Stoichiometry Indicator
Titration curves Logarithmic Measures pH & pK_a
Primary Standard

Techniques:

Titration / Back-titration Standardization
pH Measurement Weighing by Difference
Preparing solutions of Precise Concentration
in a Volumetric Flask

Apparatus:

Analytical Balance Buret Pipet
pH Meter Volumetric Flask



3

Organization of this Pre-lab Lecture

Titration	Procedure 2 – NaOH Standardization
End Points & Equivalence Points	Backtitration
Strong and Weak Acids	Calculations 2 – NaOH Standardization
Dissociation Constants / pH / pK_a	What does the End Point Look Like?
Standardization	Procedure 3 - Unknown
Stoichiometry of our Reactions	Calculations 3 - Unknown
Acid/Base Indicators	Procedure 4 – Measure pH / Calculate K_a / pK_a
How pH varies when NaOH is added to Acetic Acid	Analysis of Errors - Standardization
Procedure 1 – Preparation of KHP Solution	Analysis of Errors - Unknown
Calculations 1 – KHP Solution	Pitfalls to Avoid

5

TITRATION


A reaction conducted by slow addition of a precisely measured volume of a **reagent** solution of known concentration to

a fixed amount of another **substance** with which **the reagent reacts** until

a **SIGNAL** indicates that a significant chemical change has occurred.

We call the appearance of the signal the **END POINT** of the titration.

The **SIGNAL** is often a **COLOR CHANGE**, but may be an observable change in another property.



5

END POINTS & EQUIVALENCE POINTS

The point at which an observable **signal** occurs is called **END POINT**.


Our interest is the point at which a **stoichiometric** amount of the reagent in the buret has been consumed by the reagent in the beaker.

This is called the **EQUIVALENCE POINT**

Stoichiometric: having consumed an equivalent number of moles.

e.g., in the titration of a fixed amount of HCl with NaOH

$HCl + NaOH = NaCl + H_2O$
Equivalence point is: mol NaOH added = mol HCl



6

END POINTS & EQUIVALENCE POINTS

For a **signal** to be useful, it must be an accurate indication of completion of reaction, I.e, it must signal the **EQUIVALENCE POINT**

Signal agents (**indicators**) are chosen so that, as closely as possible,

$$\text{END POINT} = \text{EQUIVALENCE POINT}$$


7

STRONG AND WEAK ACIDS

Acids and bases can be characterized by the extent to which they **dissociate** in solution



Hydrochloric Acid is a **STRONG** acid




Acetic Acid is a **WEAK** acid



8

DISSOCIATION CONSTANTS / pH / pK_a

A **QUANTITATIVE** measure of strength or weakness of an acid (or base) is its **DISSOCIATION CONSTANT, K_a**.
For the reaction



The acid dissociation constant, **K_a**, is defined as

$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

Both [H⁺] and K_a vary over many orders of magnitude.

So, it is convenient to represent them by their logarithms instead.

$$\text{pH} = -\log_{10} [\text{H}^+] \quad [\text{H}^+] = 10^{-\text{pH}}$$

$$\text{pK}_a = -\log_{10} K_a \quad [K_a] = 10^{-\text{pK}_a}$$

9

STANDARDIZATION


The analytical determination of the purity or concentration of a substance through its reaction with a substance of **certified* composition and purity (Primary Standard)**

* NIST, US Pharmacopeia, American Chemical Society, etc.

Chemical Grades: USP, NF, Reagent, Primary Standard, Secondary Standard

This **STANDARDIZATION** involves reaction of the **base**, aqueous **NaOH**, with a **PRIMARY STANDARD**, the **weak acid**:


potassium hydrogen phthalate (KHP)



10

THE STOICHIOMETRY OF OUR REACTIONS


KHP is a monoprotic weak acid. One available proton.



The stoichiometry of the reaction is: **pK_a is 5.4**



Acetic acid, **CH₃COOH**, is also a monoprotic weak acid.




It's **pK_a is 4.7**

11

ACID/BASE INDICATORS

An Acid/Base Indicator is an organic dye whose color is different in solutions of different pH.

Indicators are often weak acids*, e.g.,



(The equilibrium lies to the left in acidic solutions.)

We use **PHENOLPHTHALEIN** which is **COLORLESS** in acidic solutions and **PINK** in basic solutions.

The color change to **pink** occurs when pH increases from < 9 to > 9.

* We use amounts of indicator sufficiently small so as not to interfere with the acid we are titrating.


12

How does the pH vary when NaOH is added to Acetic Acid

Graphs which show the dependence of pH on the volume of base added to an acid are called **Titration Curves**.

Titration curve for Acetic Acid with NaOH looks like this

13



PROCEDURE 1 - PREPARE KHP

It does not matter in which order you

- Prepare KHP & Standardize NaOH
- Titrate Unknown
- Determine pH of unknown

A. STANDARDIZATION

Prepare solution of known concentration of primary standard, **KHP**

- Weigh sample **BY DIFFERENCE***
- Dissolve fully in **ERLENMEYER**
- Bring to total volume in **VOLUMETRIC FLASK****

**** Follow manual procedure**

Must transfer ALL the solution

15

CALCULATIONS - KHP SOLUTION

Weight of vial + KHP	15.4371 g
Weight of vial + remaining KHP	12.3495 g
Weight of KHP transferred	3.0876 g
Volume of KHP Solution	0.2500 L

Follow the procedure described in the lab manual for the use of the volumetric flask!

250 mL

16

CALCULATIONS - KHP SOLUTION

Weight of vial + KHP	15.4371 g
Weight of vial + remaining KHP	12.3495 g
Weight of KHP transferred	3.0876 g
Volume of KHP Solution	0.2500 L
Molarity of KHP Solution:	
$(3.0876 \text{ g} / 204.22 \text{ g/mol}) \cdot 0.2500 \text{ L} = 0.06047 \text{ M}$	

4 Sig Figs

0.015119 Mols of KHP

17

PROCEDURE 2 - NaOH STANDARDIZATION

Determine **concentration** of stock NaOH Solution (**Nominally 0.1 M**)

Titrate

measured volumes of standard KHP solution (known concentration) with measured volumes of NaOH Solution (unknown concentration)

(Delivered from buret 1)

(Delivered from buret 2)

18

Having determined concentration of KHP solution, calculate volume of ~ 0.1 M NaOH needed to react with 40 mL of KHP.

40 mL of our KHP solution (0.06047 M) contains:
 40 mL X 0.06047 mmol / mL = 2.4 mmol of KHP
 That will require 2.4 mmol of NaOH
 What volume of 0.1 M NaOH will contain 2.4 mmol of NaOH?
 $2.4 \text{ mmol} = 0.1 \text{ (mmol/mL)} \times V \text{ (mL)}$
 $V = 2.4 / 0.1 = \sim 24 \text{ mL}$

We now know when to begin to add the NaOH from the buret slowly!

19

BACKTITRATION

In the STANDARDIZATION, BOTH REAGENTS are delivered by a BURET
 ∴ Can "BACKTITRATE"
 I.e., if end point is overshoot, can recover by adding more KHP and continuing titration, e.g.,

1. Add measured volume of KHP (~40mL)*
2. Titrate to Phenolphthalein Endpoint
3. If you overshoot endpoint, add more KHP
4. Titrate to Phenolphthalein Endpoint again

* Instead of ~35 mL label the burets.

Do at most 3 titrations for the standardization

20

CALCULATIONS - NaOH STANDARDIZATION

Molarity of KHP Solution: 0.06047 M

KHP buret reading, final	37.44 mL	
KHP buret reading, initial	3.68 mL	
Volume of KHP titrated	33.76 mL	From Part 1
NaOH buret reading, final	28.73 mL	
NaOH buret reading, initial	4.52 mL	
Volume of NaOH used	24.21 mL	
mmol of KHP titrated		
33.76 mL X 0.06047 M =	2.042 mmol	
mmol of NaOH used	2.042 mmol	Stoichiometry is 1 to 1
Molarity of NaOH		
2.042 mmol / 24.21 mL =	0.08435 M	

21

Suppose 2nd and 3rd titrations produce 0.08592 M and 0.08539 M for the NaOH concentration.

Are our three values for NaOH concentration in reasonable agreement?


Average =
 $(0.08435 + 0.08592 + 0.08539) / 3 = 0.08522 \text{ M}$

Avg Dev =
 $(0.00087 + 0.00070 + 0.00017) / 3 = 0.00058 \text{ M}$

Pct Dev =
 $100 \times 0.00058 / 0.08522 = 0.68\%$


22

WHAT DOES THE END POINT LOOK LIKE?




Before end point

KHP & Unknown solutions are colorless.




At end point

Just past end point



Way past end point



23

PROCEDURE 3 - UNKNOWN

Unknowns are solutions of acetic acid in water.
 We titrate aliquots* of (undiluted) unknown.
 In this exercise, an aliquot is a 5 mL pipet-ful (5.00 ± 0.01) of solution'

- Dilute unknown with distilled water (~ 40 mL)
- Add Phenolphthalein as indicator
- Titrate unknown to phenolphthalein end point.

Remember: you cannot backtitrate in this case

Since aliquots are identical, you can track the precision of titrations by calculating the percent error in the volumes of NaOH used.

* aliquots are fixed, repetitive fractions of a solution

24

CALCULATIONS - UNKNOWN (36)

Volume of Unknown	5.00 mL		
Concentration of NaOH solution	0.08522 M		
NaOH buret, final	22.47 mL		
NaOH buret, initial	3.15 mL		
Volume NaOH used	19.32 mL		
mmol of NaOH used	$19.32 \times 0.08522 = 1.646$ mmol		
mmol of Acetic Acid titrated	= 1.646 mmol		
Acetic Acid Concentration	$1.646 / 5.00 = 0.329$ M		

Stoichiometry is 1 to 1

3 Sig Figs

From Standardization

CALCULATIONS - UNKNOWN (36)

Volume of Unknown	5.00 mL	5.00	5.00
Concentration of NaOH solution	0.08522 M		
NaOH buret, final	22.47 mL	21.16	23.72
NaOH buret, initial	3.15 mL	2.37	4.22
Volume NaOH used	19.32 mL	18.79	19.50
mmol of NaOH used	$19.50 \times 0.08522 = 1.646$ mmol	1.601	1.662
mmol of Acetic Acid titrated	= 1.646 mmol	1.601	1.662
Acetic Acid Concentration	$1.662 / 5.00 = 0.329$ M	0.320	0.332

CALCULATIONS - UNKNOWN

Avg Conc of Acetic Acid = $(0.329 \text{ M} + 0.320 \text{ M} + 0.332 \text{ M}) / 3 = 0.327 \text{ M}$

Avg Deviation = $(0.002 + 0.007 + 0.005) / 3 = 0.005 \text{ M}$


Pct Deviation = $100 \times 0.005 / 0.327 = 1.5\%$

PROCEDURE 4 - MEASURE pH

To determine the K_a (and pK_a) of your unknown, you need to determine the pH of the unknown.

2. **Measure pH** of *undiluted* unknown to determine its H^+ concentration using a pH METER

Electronic device designed to measure hydrogen ion concentration in aqueous solutions



Since HOAc (aq) ionizes to give H^+ and OAc^- , the pH gives us the concentration of the **dissociated HOAC** in our vinegar, $[H^+] = [OAc^-]$

CALCULATE pH / pK_a

The titration gives us the **TOTAL CONCENTRATION** of HOAc in the unknown.

Suppose result of our titration is $[HOAc] = 0.124 \text{ M}$ and, the **measured pH** = 2.7

I.e., $[H^+] = 10^{-2.7} = 2.0 \times 10^{-3}$

From pH $\rightarrow [H^+]$

$[OAc^-] = [H^+]$

Total Acid $\rightarrow [HOAc]$

$K_a = \frac{[H^+][OAc^-]}{[HOAc]} = \frac{(2.0 \times 10^{-3})^2}{0.124} = 3.2 \times 10^{-5}$

$pK_a = -\log_{10}(3.2 \times 10^{-5}) = 4.5$

Should subtract 2.0×10^{-3} But to 2 sig figs, can ignore


Read & Record Burets

24.64 mL

Read & Record Weights

4.6427 g

Begin each titration with buret reading between 0.00 and 5.00 mL



ANALYSIS OF ERRORS

STANDARDIZATION - Precision

1. WEIGH ~3 g of solid KHP using the analytical balance (± 0.0004 g)
Precision: $100 \times 0.0004 / 3.0000 = 0.01 \%$

2. PREPARE A SOLUTION of KHP using 250 mL volumetric flask (± 0.05 mL)
Precision: $100 \times 0.05 / 250.00 = 0.02 \%$

3. TITRATE measured volumes of KHP and NaOH using burets (± 0.05 mL)
Precision: KHP $100 \times 0.05 / 30.00 = 0.2 \%$
NaOH $100 \times 0.05 / 30.00 = 0.2 \%$

Concentration of NaOH determined with **PRECISION** of approximately $0.01 + 0.02 + 0.2 + 0.2 = 0.43 \%$ ³¹

ANALYSIS OF ERRORS

TITRATION OF UNKNOWN - Precision

1. TRANSFER aliquot of unknown using a 5 mL transfer pipet (Intrinsic error: ± 0.01 mL)
Precision: $100 \times 0.01 / 5.00 = 0.2 \%$

2. TITRATE aliquot using standardized NaOH, again using a buret (Intrinsic error: ± 0.05 mL)
Precision: $(\text{Buret}) 100 \times 0.05 / 30.00 = 0.2 \%$

3. STANDARDIZED NaOH concentration
Precision: $= 0.4 \%$

OVERALL PRECISION of determination of concentration of unknown = $0.2 + 0.2 + 0.4 = 0.8 \%$

32

PITFALLS TO AVOID

- Incomplete transfer of solid KHP into flask (caused by weighing using a watch glass or paper)
Weigh by difference!
- Failure to dissolve KHP completely
Make sure no KHP particles remain undissolved
- Incomplete transfer of KHP solution from Erlenmeyer to volumetric flask
Transfer completely and rinse with distilled water from wash bottle
- Failing to bring KHP solution to correct volume
Use dropper for last 0.5 mL!

33

- Incomplete mixing of KHP solution in volumetric flask
Mix thoroughly both before and after bringing to final volume
- Improper use of the pipet
Make sure you use the pipet as instructed
- Contamination of unknown by allowing contact with syringe while pipeting
Discard any sample for which this happens
- Over-titrating unknown
Stop at earliest detectable pink color which persists for ≥ 30 seconds

GRADE ON THIS EXERCISE WILL DEPEND ON both Accuracy and Precision

34

How Serious is Overtitrating?

Suppose you over-titrate by 5 drops
I.e., you go past the end point by 5 drops (0.25 mL).
How large an error in accuracy results?

- Unknowns range in concentration from 0.5 M to 1.0 M.
- 5 mL of the unknown contains between 2.5 and 5.0 mmol of acetic acid.
- That will require between 25 and 50 mL of ~ 0.1 M NaOH

The error represented by 0.25 mL will range between:
 $100 \times 0.25 / 25 = 1 \%$ to $100 \times 0.25 / 50 = 0.5 \%$

Errors in identifying the end point are not a major contribution to errors in the accuracy of the concentration of the vinegar unknown.

35