

## ASPIRIN PURITY BY pH TITRATION

Last Update: 3/15/2009 1:14 PM

### Purpose:

To study pH variation during an acid/base titration and use this to determine the purity of your synthesized aspirin

### Concepts:

pH	Nernst Equation
Calibration and use of the pH meter	
Molar Mass	Effective Molar Mass
Titration	Stoichiometry
Equivalence Point	Titration of a mixture

### Techniques:

Weighing by Difference      Using a pH Meter  
Graphing

### Apparatus:

pH Meter                      Buret

**NERNST EQUATION** describes dependence of electrode voltages on concentration and temperature

For normal pH electrodes,

$$E = E^{\circ} - \frac{2.303 RT}{n_A e} \log_{10} [H^+]$$

where:  $E$  = **VOLTAGE DIFFERENCE** between reference and indicator electrodes

$E^{\circ}$  = **CONSTANT** depending only on **TYPES of ELECTRODES**

$T$  is **Absolute** Temperature

$R$  is gas constant

$n_A$  is Avogadro's number

$e$  is absolute value of electron charge

### CALIBRATION ( STANDARDIZATION )

Like all measuring devices, must calibrate pH METER before using it.

Immerse pH electrode into one or more (buffer) solutions of accurately known pH.

Set **Calibration Control** on meter so reading agrees with pH of standard buffer

Remember to rinse  
electrode after  
standardization

**SAMPLE TABLE FOR pH TITRATION**

Buret Rdg	Cum Vol	pH	Incr	Buret Rdg	Cum Vol	pH	Incr
4.35	0.00	2.94		27.10	22.75	6.35	
6.46	2.11	3.15		27.20	22.85	7.00	
8.72	4.37	3.31		27.27	22.92	8.93	1 drop
....	....	....		27.32	22.97	9.40	
24.81	20.46	4.12		27.39	23.04	10.49	
25.72	21.37	4.70		27.46	23.11	11.57	
26.39	22.04	5.20		....	....	....	....

Start using 1 drop increments about 2 mL before calculated end point.

### Phenolphthalein titration:

0.2433 g requires 21.67 mL of NaOH

How much should I weigh to use 25 mL?

$$0.2433 \text{ g} / 21.67 \text{ mL} = X / (25 \pm 3)$$

I actually weigh out 0.2742 g

**pH Titration:** 0.2742 g will require:

$$21.67 \times 0.2742 / 0.2433 = 24.42 \text{ mL of NaOH}$$

### WHAT WE ARE LOOKING FOR


Your synthesized aspirin (ASA) most likely contains some unreacted salicylic acid (SA) -  
(Did the  $FeCl_3$  test show any violet color?)

If the sample is dry, then it can contain only:

	Molar Mass	$pK_a$
Aspirin	180.15	3.50
Salicylic Acid	138.12	2.99

Titration of these acids *individually* with NaOH should allow them to be easily distinguished.

From Last Semester HOW DOES THE STRENGTH ( $pK_a$ ) OF AN ACID AFFECT ITS TITRATION CURVE?



### What about a mixture of the two?

End point in acid-base titration is characterized by rapid increase in pH

This will only occur after **both** acids are consumed, the stronger acid being consumed first.

Total amount of NaOH is the amount necessary to react with both acids.

Suppose our sample weighs  $w_{SA+ASA}$  and contains  $n_{SA}$  moles of SA and  $n_{ASA}$  moles of ASA

From the **known** number of moles of NaOH used in the titration,  $n_{NaOH}$ , we can determine

$$\Rightarrow n_{NaOH} = V_{NaOH} \times M_{NaOH} = n_{SA} + n_{ASA}$$

the sum of the moles of acid in the original sample.

But the (**known**) weight of the original sample was:

$$w_{SA+ASA} = n_{SA} \times 138.1 + n_{ASA} \times 180.1$$

$w_{SA} \qquad \qquad w_{ASA}$

We have two (linear) equations in two unknowns

These can be solved simply. The EMM formula

$$\%ASA = \frac{100 (EMM - 138.1)}{180.1 - 138.1} = \boxed{\phantom{000}}$$

is one way of representing the solution using the

$$\text{quantity, EMM} = \frac{w_{SA+ASA}}{n_{SA} + n_{ASA}} = \frac{w_{SA+ASA}}{n_{NaOH}}$$

**\* Equations involving EMM apply ONLY when EMM is between 138.1 and 180.1**

In normal (indicator only) titrations, can determine (only) number of available moles of acid in a sample - i.e., EMM

e.g., if 25.00 mL OF 0.100 M NaOH is required to react with 425 mg of a sample, the number of moles of acid in the sample is

$$25.00 \text{ mL} \times 0.100 \text{ mmol / mL} = 2.50 \text{ mmol}$$

and sample behaves like an acid of EFFECTIVE MOLAR MASS (EMM)

$EMM = 425 \text{ mg} / 2.50 \text{ mmol} = 170 \text{ mg / mmol}$   
Consistent with sample being largely ASA (MM = 180) with some SA impurity (MM = 138)

**BE SURE TO****LEAVE ELECTRODE IN SOLUTION WHEN DOING pH TITRATION**

Improves reliability and saves time

**STIR THE CONTAINER** into which you are titrating to obtain good mixing.

Otherwise, you measure pH of only a small, local part of the solution

**RINSE THE pH ELECTRODE WITH DISTILLED WATER and DRY IT CAREFULLY** each time you immerse it into a **NEW** liquid

Otherwise, pH electrode can still "see" the last solution with which it was in contact

**Make abscissa of graph CUMULATIVE VOLUME of added NaOH**

Remember to subtract the initial buret reading

**Record NaOH STOCK SOLUTION CONCENTRATION**

You need it to calculate EMM's

**SAVE REMAINDER OF SYNTHESIZED ASPIRIN FOR FUTURE EXERCISE****pH METERS MAY NEED OCCASIONAL RECALIBRATION - between users****SUMMARY OF PROCEDURE**

- 1.) Working in **ASSIGNED PAIRS**, perform **NORMAL TITRATIONS** on samples of:  
**AUTHENTIC SALICYLIC ACID (SA)**  
**AUTHENTIC ACETYSALICYLIC ACID (ASA)**  
You will be given vials of pure SA and ASA
- 2.) Compute EMM's for both of the above
- 3.) Working **ALONE** from here on, titrate **YOUR SAMPLE OF ASA** with standardized NaOH, using phenolphthalein as an indicator

Agreement of your EMM's with expected values for SA and ASA is a measure of your titration technique.

- 4.) From this, determine EMM **and** the **AMOUNT** of **YOUR SAMPLE OF ASA** needed to consume 25 mL of the same NaOH solution
- 5.) Weigh out your ASA sample and calculate how much NaOH it will require.
- 6.) Set up pH meter, check calibration
- 7.) **TITRATE YOUR SAMPLE OF ASA** with standardized NaOH, using **pH METER**

Weigh all samples on **ANALYTICAL BALANCE** by **DIFFERENCE**

**DO NOT** TRANSFER ANY COMPOUND TO BALANCE PAN!

Even though you share pH meters,

**Don't do simultaneous pH titrations**

Datasheet and graph for this exercise may be turned in **at the beginning of the next lab**

Do NOT record pH titration data on data sheet - only in lab notebook

**Calculations****Regular titration of synthesized ASA**

Weight of sample	205.6 mg
Volume of NaOH	22.61 mL
Conc of NaOH	mmol NaOH 0.05246
EMM =	

$$205.6 / (22.61 * 0.05246) = 177.4 \text{ mg}$$

**What weight is required for 25 ± 3 mL of NaOH?**

$$(25 \pm 3) * 205.6 / 22.61 = (25 \pm 3) * 9.093 = 227 \pm 27 \text{ mg}$$

Suppose you **actually** weigh 211.5 mg of sample.

**Where will the end point be?**

$$22.61 * 211.5 / 205.6 = 23.26 \text{ mL NaOH}$$

Should begin 1 drop increments at ~21 mL

### Calculations - Continued

After doing the pH titration on your sample, you will have three quantities that let you determine the purity of your synthesized ASA:

Your  $EMM_{sa}$  for SA and

Your  $EMM_{asa}$  for ASA (from the regular titrations),

Your  $EMM_s$  for your synthesized aspirin.


You should calculate your purity using these numbers, i.e.,

$$\%ASA = \frac{100 (EMM_s - EMM_{sa})}{EMM_{asa} - EMM_{sa}}$$


**YOUR  
VALUES**

### REMINDERS ABOUT GRAPHING (SUPL-004)

- Use rational number of boxes per unit (1, 2, 5, 10, 20, ..... ) for pH and volume
- Arrange graph so that maximum total area of graph paper is utilized
- Draw a smooth curve through the experimental points
- Label equivalence (end) point and half-titration point clearly
- Interpolate values of these points with precision (significant figures) consistent with your plot
  - *You should be able to read the ordinate and abscissa to at least the nearest 0.1 mL or 0.1 pH unit*



You must plot the pH vs volume of NaOH added **BY HAND** -  
not on a computer



pH axis does not need to start at 0

### NEXT WEEK

Next Week's Exercise uses same pH TITRATION TECHNIQUE to IDENTIFY AN

UNKNOWN ACID

*TEST EXERCISE*

READ SUSB-014  
DO PRELAB