



**OPTIONAL:** Weigh a second sample of the computed amount of your acid and repeat steps 3 - 6.

Do this only after everyone has finished one pH titration.

*May wish to repeat normal titration*

You obviously must work alone, but pH meters will be shared by as many as three students

Titration curve permits the determination of two characteristics of the unknown acid


1) END POINT determines the EFFECTIVE MOLAR MASS (EMM)

2) pH at HALF EQUIVALENCE POINT determines pK<sub>a</sub>

**YOU MUST NOT SHARE pH METERS DURING A TITRATION**

**THE SECOND USER MUST WAIT UNTIL THE FIRST USER HAS FINISHED THE pH TITRATION**

### The Half Titration Point




At the

$$[H^+] = K_a \frac{[HA]}{[A^-]}$$

In the form:

$$pH = pK_a + \log \frac{[A^-]}{[HA]}$$


This is the famous *Henderson-Hasselbalch equation*.



### EMM (from end point):

E.g., If 188.6 mg of your unknown requires 26.37 mL of 0.0533 M NaOH to reach the end point:

Titration curve - 188.6 mg of Weak Acid with 0.0533 M NaOH



### EMM (from end point):

E.g., If 188.6 mg of your unknown requires 26.37 mL of 0.0533 M NaOH to reach the end point:

$26.37 \text{ mL} \times 0.0533 \text{ mmol/mL} = 1.41 \text{ mmol}$  of NaOH were used, and

$$\text{EMM} = 188.6 \text{ mg} / 1.41 \text{ mmol} = 134 \text{ mg/mmol}$$

pK<sub>a</sub>: from pH (at half equiv pt)  $\approx$  pK<sub>a</sub>

In this example, half equivalence point was at:

$$26.37 \text{ mL} / 2 = 13.19 \text{ mL}$$

pH at that volume of added NaOH was 2.2,

pK<sub>a</sub> of acid is:

$$pK_a = 2.2$$

The titration curve has produced:

$$\text{MOLAR MASS} = 134 \text{ g / mol}$$

$$\&$$

$$\text{p}K_a = 2.2$$

To what acid do these data correspond?

What are the errors in these two numbers and what contributes to the errors?

TABLE I: WEAK ACIDS

Glutaric acid	$\text{H}_2\text{C}_5\text{H}_6\text{O}_4$	132.1	4.3
dl-Malic acid	$\text{HC}_4\text{H}_4\text{O}_5$	134.1	3.4, 5.0
Potassium dihydrogen phosphate	$\text{KH}_2\text{PO}_4$	136.1	7.2, 12.7
Potassium bisulfate	$\text{KHSO}_4$	136.2	1.9
Sodium bisulfate hydrate	$\text{NaHSO}_4 \cdot \text{H}_2\text{O}$	137.4	1.9
Sodium dihydrogen phosphate	$\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$	138.0	7.2, 12.7

Our numbers are 134 2.2

### WHAT CAN CAUSE ERRORS IN THIS EXERCISE?

#### 1. SAMPLE WEIGHT (~200 mg)

Weigh **BY DIFFERENCE** !!!!!!!!!!!!!

Suppose you weigh on watch glass and lose 10 mg of sample in transfer

$$-10 / 200 = -5\% \text{ ERROR}$$

#### 2. TITRATION (~25 mL)

Miss End Point by  $\pm 0.50$  ml (10 drops)

$$\pm 0.50 / 25 = \pm 2\% \text{ ERROR IN EMM}$$

$$\pm 0.25 / 12.5 = \pm 2\% \text{ ERROR IN HALF TITRATION VOLUME}$$

&

$$\pm 1\% \text{ ERROR IN p}K_a$$



While primary purpose of part 1 (titer determination) is to determine appropriate weight of sample to use for the pH titration,

first titration *also* provides a good estimate of acid EMM. Can use it to verify EMM obtained from pH titration.

If EMM's from regular titration and pH titration differ by 10% or more, it is likely that you have not plotted pH vs the appropriate values for the delivered volume.

**YOU ARE URGED TO USE THE SAME BALANCE FOR ALL MASS MEASUREMENTS**

### The Half Equivalence Point

How good is our approximation that **pH at half-equivalence point = pK<sub>a</sub>**?

For the concentrations used in this exercise and the pK<sub>a</sub>'s of the unknown acids, the error in the approximation

*pH = pK<sub>a</sub> at 1/2 of the volume of base delivered at the equivalence point*

**is less than 10%**

For pK<sub>a</sub>'s > 4, it is less than 1%

A web page deals extensively with this approximation

### SAMPLE DATA SHEET

Mass of vial + Acid	14.6705 g
Mass of vial - sample	14.5581 g
Mass of Acid Sample (g)	0.1124 g
Mass of Acid Sample (mg)	112.4 mg FROM CONTAINER
Vol of NaOH @ Equiv Pt	24.5 mL FROM GRAPH
Concentration of NaOH	0.0533 M
mmol of NaOH	1.31 mmol FROM GRAPH
mmol of Acid	1.31 mmol
Molar Mass of Acid	85.83 mg/mmol
Vol NaOH @ Half Equiv	12.3 mL
pH @ Half Equiv	5.0
pK <sub>a</sub> of Acid	5.0

The Data are Consistent with Crotonic Acid (MM = 86.1; pK<sub>a</sub> = 4.70) [TABLE]

**NEXT EXERCISE**

**IDENTIFICATION OF PLASTICS**

READ **SUSB-013**

And

**DO** PRE-LAB from **SUSB-013**